


### Flow Injection Ammonium Analyzer

The automated analyser for measuring low concentrations of ammonium in water. The FloRRia can hold as much as 100 sample tubes at a time performing the analyses unattended.



Only a NaOH-solution and demi-water is needed to perform analyses. Sample rate is approx. 12 samples per hour. Concentrations of 5 ppb (minimum) to 30.000 ppb (maximum) may be measured.



Interfacing through RS-232 to a host computer running easy-to-use dedicated software enabling calibration, operation, setup, control and archiving.

#### Operation principle

Aspirating from a 0.1 ml tube fully automated, the sample is mixed with NaOH and subsequently passes the membrane inside the detector block. Part of the gas-phase ammonia diffusion passes through the membrane and is there added to demi-water. The conductivity of the solution is compared with the conductivity of neutral demi-water with that defining the ammonia concentration in the sample.



#### Summary of specifications

<i>Measuring characteristics</i>	
Dynamic range	5 ppb – 30,000 ppb
Resolution	5 ppb
Accuracy	5 %
<i>Sampling characteristics</i>	
Throughput	12 samples/hour
Sample tubes	100 max
Tube size	Ø 12 - 13.5 mm L=70 - 80 mm
Sample volume	1.5 ml
<i>Environmental requirements</i>	
Temperature	15 – 30 °C
Relative humidity	10 – 90 %
<i>Power</i>	
Power consumption	12 VDC, 36 Watt max.
Power converter input	100-250 VAC; 47/63 Hz; 0.8A
output	12 VDC; 3 A
<i>Dimensions and weight</i>	
Overall dimensions	620 x 270 x 400mm (LxWxH)
Weight	9 kg (excl. reagents)

All specifications subject to change without prior notice

[www.mechatronics.nl](http://www.mechatronics.nl)